



ALTERNATING PRESSURE PUMPS

- COMPLETE WITH MATTRESS
- HELPS PREVENT PRESSURE SORES
- DOUBLE-ISOLATED MOUNTING FOR QUIETEST OPERATION



A



B

A Model #2600

**B Model #2650 (Adjustable)
(Model #2675 - Replacement Mattress)**

These state-of-the-art alternating pressure pad and pump systems change pressure points on a regular cyclical basis, thereby bringing relief to large capillary areas and permitting more beneficial circulation. They help prevent pressure sores and relieve unnecessary suffering. Their double-isolated mounting ensures the quietest operation. All connections and controls are convenient and simple to use. The pad system contains two sets of air cells which alternately inflate and deflate every two and one half minutes. To minimize movement, full inflation of one set of cells occurs before the large alternate set is deflated. Large cell interlocking support patterns permit a high volume of air exchange and maximum vertical excursion between relief and support modules. Electronically welded heavy gauge premium vinyl, which is reusable, resists tearing and accidental punctures.

Care & Use: Do not use in the presence of smoking materials or open flame (air flowing through pad will support combustion).

Specifications: Pad: 100% welded heavy gauge vinyl.

	Model#2600	Model#2650
Pressure Range	2.1 PSI	1.1 to 2.1 PSI
Air Output	110V: 4.5 liters / 220V: 5.0 liters	110V: 4 to 5 liters / 220V: 5 to 6 liters
Weight	3.4 lbs	3.5 lbs
Power	115V/60Hz / 220V/50Hz / 100V/50Hz, 60Hz	
Cycle Time	110V: 2.5 minutes / 220V: 3.0 minutes	
Pump Dimensions	9 3/4" x 5 1/8" x 3 7/8"	
Mattress Dimensions	80" x 36"	
Power Cord	10'	
CSA / UL Approval	✓	✓

